



Mouse/Rabbit Probe HRP Labeling Kit with DAB Brown

Product Name: Mouse/Rabbit Probe HRP Labeling Kit with DAB Brown

Catalog No.: TAHC03D

Form:

| Catalog No. | Size |
|-------------|-------|
| TAHC03D-15 | 15ml |
| TAHC03D-50 | 50ml |
| TAHC03D-100 | 100ml |

INTENDED USE: For Research Use Only. Not for use in diagnostic procedures.

DESCRIPTION: The reagents in this kit constitute a biotin-free immunoenzymatic antigen detection system. This technique involves the sequential incubation of the specimen with an unconjugated primary antibody specific to the target antigen, a secondary antibody- HRP conjugate which reacts with mouse and rabbit primary antibody, and substrate-chromogen (DAB).

APPLICATIONS: Immunohistochemistry (IHC)

| KIT CONTENTS: | Description | Format | Recommend time |
|---------------------------------|---|--------------|--|
| Hydrogen Peroxide Block | 3% hydrogen peroxide solution with less than 0.1% sodium azide | Ready to Use | 10-20 minutes. |
| Immunoblock | PBS solution, pH 7.6, with 0.5% BSA, and less than 0.1% sodium azide. | Ready to Use | 10-30 minutes. |
| Mouse/Rabbit Probe HRP Labeling | Polymer conjugated HRP Reactive with rabbit primary antibodies. | Ready to Use | 30-45 minutes. |
| DAB Chromogen (20x) | 3,3' Diaminobenzidine (DAB) chromogen | Concentrate | Diluted with DAB Buffer (chromogen 1part : buffer 19 parts) |
| DAB Buffer | Substrate buffer, pH 7.5 with hydrogen peroxide | Ready to Use | 1-10 minutes. |
| Hematoxylin | Hematoxylin | Ready to Use | 0.5-1 minutes. |

STAINING PROTOCOL:

1. Deparaffinize and rehydrate formalin-fixed paraffin-embedded tissue section.
2. Add enough drops of Hydrogen Peroxide Block to cover the sections. Incubate for 10 minutes. Wash 2 times in PBS buffer.
3. Perform appropriate pretreatment if required. Wash 3 times in PBS buffer.
4. Apply Immunoblock and incubate for 10-30 minutes at room temperature to block



- non-specific background staining. Wash 2 times in PBS buffer.
5. Apply mouse or rabbit primary antibody and incubate according to manufacturer's protocol.
Wash 3 times in PBS buffer.
 6. Apply Mouse/Rabbit Probe HRP Labeling and incubate for 30 minutes at room temperature. Wash 3 times in PBS buffer.
 7. Add 50 μ l DAB Chromogen to 950 μ l of DAB Buffer (20X dilution), mix by swirling and apply to tissue. Incubate for 1-10 minutes. Rinse 4 times in PBS buffer.
 8. Rinse slide in tap water.
 9. Apply Hematoxylin counterstain according to manufacturer's instructions (optional).
 10. Dehydrate if required and cover slips.

STORAGE & STABILITY:

Store at 2-8°C. Do not freeze. The user must validate any other storage conditions. The reagents must be returned to the storage conditions identified above immediately after use. When properly stored, the reagent is stable to the date indicated on the label. Do not use the reagent beyond the expiration date. There are no definitive signs to indicate instability of this product; therefore, positive and negative controls should be tested simultaneously with unknown specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact Technical Support at info@biotna.net or your distributor service.