



Aqua mount

Product name: Aqua mount

Catalog No. : TA00M2

Form:

Catalog No.	Size
TA00M2 - 50	50 mL
TA00M2 - 100	100 mL
TA00M2 - 200	200 mL
TA00M2 - 500	500 mL

INTRODUCTION

Aqua Mount is a non-permanent aqueous mounting medium formulated for coverslipping directly from aqueous solutions. It is non-fluorescing and has an antifade component to increase the viewing time of the specimen. Use Aqua Mount with most fluorescent dyes and stains including DAB, Alkaline Phosphatase Fast Red, AEC (aminoethylcarbazole) and a variety of other chromogens to enhance and retain fluorescent intensity.

Aqua Mount can be used for frozen sections, fat stains, chromogens for immunohistochemistry and in situ hybridization as well as other applications requiring a water soluble mounting medium. Also, useful in avoiding leaching or de-staining of enzymatic staining protocols like TRAP. Using Aqua Mount, the coverslip will dry in 24 to 48 hours depending on the amount of medium used. Aqua Mount solidifies under coverslip on the microscope slide after 24 hours with the following laboratory conditions and ambient temperature; relative humidity in winter time at greater than 20% and relative humidity in the summer time of up to 50%. Painting the edges of the coverslip is unnecessary. When the medium dries it will form a seal.

PROCEDURE

Prepare slides as required. Prior to coverslipping rinse the slides in distilled or deionized water. The excess fluid can be drained, but blotting is not recommended. Aqua Mount is supplied in convenient 20ml dropper bottles for manual staining and in larger volumes for automated coverslippers.

Aqua Mount can be used on automated coverslippers in place of traditional solvent based coverslipping medium. The viscosity of the medium may require adjusting for proper flow rates. Allow the liquid to come to room temperature before making any adjustments for flow rate as the viscosity will change. The amount of medium dispensed will effect the drying time. Thicker sections will require more coverage with medium to assure the coverslip remains in place against the slide surface. contact Technical: info@biotna.net