



**BIOTnA**  
molecular pathology


# CATALOGUE

# Our Passion



At BIOTnA our passion is providing biomedical laboratories with the tools to improve the diagnosis, prognosis and therapies that benefit patients worldwide.

BIOTnA performs R&D, production, distribution and marketing of unique products for Immunohistochemistry (IHC), Immunocytochemistry (ICC), Fluorescent in situ hybridization (FISH), Chromogenic in situ hybridization (CISH), In Situ PCR and Polymerase Chain Reaction (PCR) technologies that meet the highest international standards for applications in Molecular Pathology, Cancer Research, Microbiology, Immunology and Genetics.



# Table of Contents

## 01 > Immunohistochemistry

---

### IHC

Mouse Probe HRP Labeling Kit.....	02
Rabbit Probe HRP Labeling Kit.....	03
Mouse/ Rabbit Probe HRP Labeling Kit.....	04
Double stain Kit.....	05
Substrate Chromogen DAB AEC Green Blue Alk-Red Alk-Blue.....	05

### CISH · FISH

BioSpot ISH Kit For DIG probe.....	07
BioSpot ISH Kit For Biotin probe.....	07
BioSpot ISH Kit with FAM For DIG Probe.....	08
BioTnA In Situ PCR kit with DAB Brown for DIG labeled probe.....	08 · 09
BioTnA In Situ RT-PCR Kit with DAB Brown for DIG labeled probe....	10 · 11

### TUNEL

TUNEL Apoptosis Assay Kit.....	12
TUNEL Apoptosis Assay Kit (FAM).....	12

## 02 > Immuno-Fluorescence


---

Fluorescence multiple stain kit.....	13
Fluorescence Blocking Reagent.....	13
Goat anti-Mouse IgG (H+L)-488.....	13
Goat anti-Mouse IgG (H+L)-594.....	13
Goat anti-Rabbit IgG (H+L)-488.....	13
Goat anti-Rabbit IgG (H+L)-594.....	13
Donkey anti-Goat IgG (H+L)-488.....	13
Donkey anti-Goat IgG (H+L)488.....	13
Goat anti-Rabbit IgG-iFluor-680.....	13

## 03 > Special stain

---





Masson Trichrome.....	14
Periodic Acid-Schiff(PAS).....	14
Alcian blue(pH2.5 · pH1.0).....	14
Nissl stain.....	14
Von Kossa.....	15
Oil red O.....	15
Toluidine blue.....	15
PicroSirius red.....	15
Alizarin Red S.....	15
Safranin O.....	16
Verhoeff-Van Gieson stain(VVG).....	16
Congo red.....	16
Prussian blue.....	16
Luxol-Fast blue.....	16



# Table of Contents

## 04 > Ancillaries for Immunohistochemistry

Hematoxylin.....	17
Nuclear blue for IHC (modified hematoxylin).....	17
Eosin.....	17
Hydrogen Peroxidase Block.....	17
Background Reducing buffer(general).....	17
Background Reducing buffer(goat).....	17
Davidson' s Fixative solution.....	17
20X Citrate Buffer, pH 6.0.....	17
20X Tris EDTA Buffer, pH9.0.....	17
Enzyme Digester I (Pepsin).....	17
Enzyme Digester II(proteinase K).....	17
Enzyme Digester III(Trypsin).....	17
Permont Mount.....	17
Aqua Mount.....	17
1X PBS, (powder).....	17
TBST.....	18
20X SSC Buffer.....	18
Antibody Diluent.....	18
Nuclear Fast Red.....	18
Fast green.....	18
Methyl Green.....	18
Goat probe.....	18
Rat probe.....	18
BioSpot ISH Kit For DIG probe.....	18
BioSpot ISH Kit For Biotin probe.....	18
Goat anti-Mouse IgG (H+L)-HRP.....	18
Goat anti-Rabbit IgG (H+L)-HRP.....	18
Donkey anti-Goat IgG (H+L)-HRP.....	18
Normal Mouse IgG.....	18
Normal Rabbit IgG.....	18
Normal Goat IgG.....	18
Hydrogen PeroxidaseBlock.....	18
Fluorescence Blocking Reagent.....	18
Background Reducing buffer.....	18
Background Reducing buffer plus.....	18



## 05 > Molecular biotechnology

---

TnA Taq DNA polymerase.....	19
Hot start Taq DNA polymerase.....	20
Hi Fidelity polymerase.....	21
Pfu DNA polymerase.....	22
2x PCR mix.....	23
2x Redy mix.....	23
2X Hi Fi mix.....	24
2X Hot start mix.....	24
2X Pfu mix.....	25
5X GC rich buffer.....	25
2x PCR mix.....	26
Probe 2X qPCR Master mix.....	27
M-MLV reverse transcriptase H- .....	28
2X One-tube RT-PCR mix.....	29
2X RT Master Mix.....	30
RNAsave solution.....	31
RNase Eraser.....	31
5X RN A Protector.....	31
100 bp DNA Ladder.....	32
100 Kb DNA Ladder.....	32
T4 DNA ligase.....	33
dNTP mix.....	34
dNTP set.....	34
RNA reagent.....	35
Genomic DNA reagent.....	35
Blue Stain.....	36
Bradford reagent.....	36
cDNA Synthesis Kit.....	37
2X qPCR Master mix.....	38
Probe 2x qRT-PCR Master mix.....	39
HyFect™ DNA Transfection Reagent....	40
DNA Safe Stains.....	41

Immunohistochemistry or Immunocytochemistry is the name of the method for localizing specific antigens in tissue or cells using antibodies, enzyme conjugates and substrate-chromogens. The antigen-antibody reaction can be evidenced with an optical microscope. The antibodies can be polyclonal or monoclonal (usually from mouse or rabbit) in origin. These non-infectious biological reagents are used by research and clinical laboratories to detect the presence of infectious agents, cancer and other proteins of interest.

BIOTnA designs several polymer detection kits for signal and double stain. We also design a special blocker as a protein blocker that can reduce the non-specific background caused by endogenous IgG and IgM avoids the mouse on mouse problems. We offer HRP been developed using a proprietary tandem hyperlabelling technology used to directly label immunoglobulins with enzymes. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranous antigens, in different types of tissues or cells.

We supply our customers with a wide array of high sensitivity and easy-to-use Substrate-chromogen Systems for HRP (DAB, AEC, HRP Green). Our Ancillary products for IHC include different Retriever Solutions for HIER, Enzymes for Tissue or Cell digestion, Diluents and Blockers, Reagents Controls, Buffers, Mounting Media and Counterstains.



### ● Mouse Probe HRP Labeling Kit



Storage condition : Store at 2-8 °C

Expired date : 18-24 months

#### Description

The Mouse Probe HRP Labeling Kit is a non-biotin, 1-step detection system that allows for the demonstration of antigens in paraffin-embedded and frozen tissue, cryostat sections, blood smears, and cytosmears. The kits have been developed using a hyperlabelling technology used to directly label immunoglobulins with enzymes. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranes antigens, in different types of tissues.

The Mouse Probe HRP Labeling Kit is suitable for use with mouse IgG and IgM primary antibodies. The Mouse Probe HRP Labeling kits are optimized for use with most primary antibodies; however, they are universal kits and therefore work equally well with pre-diluted and concentrated antibodies from different vendors' primaries.

#### Mouse Probe HRP Labeling Kit with DAB brown

Cat No	Pack size	Tests
TAHC01D-15	15 mL	150T
TAHC01D-50	50 mL	500T
TAHC01D-100	100 mL	1000T

#### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Mouse probe HRP conjugated
- DAB buffer
- DAB chromogen (20X)
- Hematoxylin

The Mouse Probe HRP Labeling Kit with DAB brown is displayed brown color. It could be dehydrated and maintained for long time.


#### Mouse Probe HRP Labeling Kit with AEC

Cat No	Pack size	Tests
TAHC01A-15	15 mL	150T
TAHC01A-50	50 mL	500T
TAHC01A-100	100 mL	1000T

#### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Mouse probe HRP conjugated
- AEC chromogen
- Hematoxylin

The Mouse Probe HRP Labeling Kit with AEC is displayed red color. It could not be dehydrated.



## ● Rabbit Probe HRP Labeling Kit



Storage condition : Store at 2-8 °C

Expired date : 18-24 months

### Description

The Rabbit Probe HRP Labeling Kit is a non-biotin, 1-step detection system that allows for the demonstration of antigens in paraffin-embedded and frozen tissue, cryostat sections, blood smears, and cytospins. The kits have been developed using a hyperlabelling technology used to directly label immunoglobulins with enzymes. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranes antigens, in different types of tissues.

The Rabbit Probe HRP Labeling Kit is suitable for use with rabbit IgG and IgM primary antibodies. The Rabbit Probe HRP Labeling kits are optimized for use with most primary antibodies; however, they are universal kits and therefore work equally well with pre-diluted and concentrated antibodies from different vendors' primaries.

The Rabbit Probe HRP Labeling Kit with DAB brown is displayed brown color. It could be dehydrated in alcohol.

### Rabbit Probe HRP Labeling Kit with DAB brown

### Rabbit Probe HRP Labeling Kit with AEC

Cat No	Pack size	Tests
TAHC02D-15	15 mL	150T
TAHC02D-50	50 mL	500T
TAHC02D-100	100 mL	1000T

Cat No	Pack size	Tests
TAHC02A-15	15 mL	150T
TAHC02A-50	50 mL	500T
TAHC02A-100	100 mL	1000T

### Contents :


- Hydrogen peroxidase block
- High effect immunoblock
- Rabbit probe HRP conjugated
- DAB buffer
- DAB chromogen (20X)
- Hematoxylin

### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Rabbit probe HRP conjugated
- AEC chromogen
- Hematoxylin

The Rabbit Probe HRP Labeling Kit with DAB brown is displayed brown color. It could be dehydrated in alcohol.

The Rabbit Probe HRP Labeling Kit with AEC is displayed red color. It could not be dehydrated in alcohol.





### ● Mouse/Rabbit Probe HRP Labeling Kit



Storage condition : Store at 2-8 °C

Expired date : 18-24 months

#### Description

The Mouse/Rabbit Probe HRP Labeling Kit is a non-biotin, 1-step detection system that allows for the demonstration of antigens in paraffin-embedded and frozen tissue, cryostat sections, blood smears, and cytosmears. The kits have been developed using a hyperlabelling technology used to directly label immunoglobulins with enzymes. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranes antigens, in different types of tissues.

The Mouse/Rabbit Probe HRP Labeling Kit is suitable for use with mouse and rabbit IgG and IgM primary antibodies. The Mouse/Rabbit Probe HRP Labeling kits are optimized for use with most primary antibodies; however, they are universal kits and therefore work equally well with pre-diluted and concentrated antibodies from different vendors' primaries.

#### Mouse/Rabbit Probe HRP Labeling Kit with DAB brown

Cat No	Pack size	Tests
TAHC03D-15	15 mL	150T
TAHC03D-50	50 mL	500T
TAHC03D-100	100 mL	1000T

#### Mouse/Rabbit Probe HRP Labeling Kit with AEC

Cat No	Pack size	Tests
TAHC03A-15	15 mL	150T
TAHC03A-50	50 mL	500T
TAHC03A-100	100 mL	1000T

#### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Mouse and Rabbit probe HRP conjugated
- DAB buffer
- DAB chromogen (20X)
- Hematoxylin

#### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Mouse and Rabbit probe HRP conjugated
- AEC chromogen
- Hematoxylin

The Mouse/Rabbit Probe HRP Labeling Kit with DAB brown is displayed brown color. It could be dehydrated in alcohol.


The Mouse/Rabbit Probe HRP Labeling Kit with AEC is displayed red color. It could not be dehydrated in alcohol.

## ● Double stain Kit

### Description

The double stain kits apply IHC technology and detect more than one target protein on the same section slide. BIOTnA designed several double stain kit in HRP system those are suitable for use with mouse and rabbit IgG and IgM primary antibodies. In our system, the powerful immunoblock could reduce the background and non-specific binding from endogenous contamination and the first primary reaction once time. The Systems included DAB brown and HRP Green chromogen. They will well define clearly two biomarkers as present in the same or different area.

- TADS01 Mouse × Rabbit
- TADS02 Goat × Mouse/Rabbit
- TADS03 Mouse/Rabbit × Mouse/Rabbit
- TADS04 Rat × Mouse/Rabbit




**Storage condition :** Store at 2-8 °C

**Expired date :** 18-24 months

Cat No	Pack size	Tests
TADS01	15 mL	150T
TADS02	15 mL	150T
TADS03	15 mL	150T
TADS04	15 mL	150T

### Contents :

- Hydrogen peroxidase block
- High effect immunoblock
- Rat or Goat probe HRP conjugated
- DAB buffer
- DAB chromogen (20X)
- HRP Green buffer
- HRP Green chromogen



## ● Substrate Chromogen

Product	Cat No	Pack size	Tests
HRP Green Chromogen Kit	TA01HG-15	15 mL	150T
HRP Green Chromogen Kit	TA01HG-50	50 mL	500T
HRP Green Chromogen Kit	TA01HG-100	100 mL	1000T
DAB Brown Chromogen Kit	TA01DB-15	15 mL	150T
DAB Brown Chromogen Kit	TA01DB-50	50 mL	500T
DAB Brown Chromogen Kit	TA01DB-100	100 mL	1000T
AEC Solution (ready-to-use)	TA01AE-15	15 mL	150T
AEC Solution (ready-to-use)	TA01AE-50	50 mL	500T
AEC Solution (ready-to-use)	TA01AE-100	100 mL	1000T

**Storage condition :** Store at 2-8 °C

**Expired date :** 18-24 months

### Contents :

- Chromogen (concentrated)
- Buffer

### INTRODUCTION

In Situ Hybridization (ISH) is a powerful technique for localizing specific nucleic acid targets within fixed tissues and cells, allowing you to obtain temporal and spatial information about gene expression and genetic loci.

### BIOTnA

ISH detection kits are designed by Chromogenic In Situ Hybridization (CISH) in formalin-fixed, paraffin-embedded tissue sections, cell samples, blood or bone marrow smears, and metaphase chromosome spreads.

### ADVANTAGES OF CISH OVER FISH

- Quick and easy interpretation of results comparable to IHC
- Simultaneous observation of tissue morphology and CISH signals
- Storage of slides at room temperature – CISH signals are permanent
- No costly fluorescent microscope needed.

#### Pretreatment Step

1. De-wax, and rehydration
2. Deperoxidase
3. Heat treatment
4. Enzyme treatment
5. Probe Hybridization 37°C, 16 hr

#### Detection Step

1. Non-specific blocking
2. Anti-DIG
3. Polymer HRP conjugated
4. DAB development
5. Counterstain 20~60s
6. Covering and reading

### ● BioSpot ISH Kit For DIG probe

#### Description

The BioSpot ISH Kit For DIG probe is designed to be used for the detection of digoxigenin (DIG) labeled probes in either formalin-fixed, paraffin-embedded tissue or cell samples by chromogenic in situ hybridization (CISH).

The kits have been developed using a hyperlabelling technology used to directly label probe with HRPs. The kit included complete pretreatment, labeling, blocking, washing and detection system. It also could be used for animal and plant samples in research filed.

TASH01D with DAB chromogen  
TASH01A with AEC chromogen

**Storage condition :** Store at 2-8 °C

**Expired date :** 18-24 months

Cat No	Pack size	Tests
TASH01D-50T	6mL	50T
TASH01D-100T	12 mL	100T
TASH01A-50T	6mL	50T
TASH01A-100T	12 mL	100T

#### **Contents :**

Hydrogen peroxidase block  
High effect immunoblock  
20X SSC buffer  
Anti-DIG labeling  
DAB or AEC chromogen  
Detection system

### ● BioSpot ISH Kit For Biotin probe

#### Description

The BioSpot ISH Kit For Biotin probe is designed to be used for the detection of Biotin labeled probes in either formalin-fixed, paraffin-embedded tissue or cell samples by chromogenic in situ hybridization (CISH).

The kits have been developed using a hyperlabelling technology used to directly label probe with HRPs. The kit included complete pretreatment, labeling, blocking, washing and detection system. It also could be used for animal and plant samples in research filed.

TASH02D with Biotin chromogen  
TASH02A with Biotin chromogen

**Storage condition :** Store at 2-8 °C

**Expired date :** 18-24 months

Cat No	Pack size	Tests
TASH02D-50T	6mL	50T
TASH02D-100T	12 mL	100T
TASH02A-50T	6mL	50T
TASH02A-100T	12 mL	100T

#### **Contents :**

Hydrogen peroxidase block  
20X SSC buffer  
DAB or AEC chromogen  
Biotin Detection system

### ● BioSpot ISH Kit with FAM For DIG Probe

#### Description

The Biospot kit F488 is designed to be used for the detection of digoxigenin (DIG)-labeled probes in either formalin-fixed, paraffin-embedded tissue or cell samples by fluorescence in situ hybridization (FISH).

Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

Cat No	Pack size
TASH01F-50	50 reactions
TASH01F-100	100 reactions

#### Contents :

##### Pretreatment kit

Heat Pre-treatment Solution, 20x

Enzyme Solution

SSC buffer, 20x

Hu DNA (+) control probe 5'-DIG

Ms DNA (+) control probe 5'-DIG

RNA (+) control probe 5'-DIG

##### Detection system

Hi-Effect Block reagent

Anti-DIG antibody

F488 labeling

DAPI Solution

### ● BioTnA In Situ PCR kit with DAB Brown for DIG labeled probe

#### Description

The BioTnA In Situ PCR Kit with DAB is designed to be used for the detection of digoxigenin (DIG)-labeled probes hybridized specific DNA sequence through PCR application in either formalin-fixed, paraffin-embedded tissue or cell samples by chromogenic in situ hybridization (CISH).

Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

Cat No	Pack size
TAISP01D-50	50 reactions
TAISP01D-100	100 reactions

#### Contents :

##### Pre-PCR treatment

Heat Pre-treatment Solution, 20x

Hydrogen Peroxidase Block

Lysis buffer

Enzyme Solution

##### PCR reagent

PCR reaction mixture

Control primer set

##### Post-PCR treatment

Fixation solution

SSC buffer, 20x

Wash buffer, 50X

##### Detection kit

Hybridization buffer

Sealing gel

Universal probe DIG labeled

Control probe DIG labeled

Hi-Effect Block reagent

Anti-DIG antibody

Polymer HRP labeling

DAB buffer

DAB chromogen, 20x

Nuclear Blue Solution

##### Slide

Negative control slide\*2

Control slide\*4


### I. primer design

### II. Sample preparing

- 1.dewax
- 2.rehydration



### III. pretreatment

- 1.De-peroxidase
- 2.Heat per-treatment
- 3.Lysis per-treatment
- 4.Enzyme pre-treatment



### IV. PCR Step

seal the cover slip in sealing gel



### V. Post-PCR Fixation

### VI. In Situ Hybridization

### VII. Detection

- 1.anti-Dig
- 2.HRP Labeling
- 3.DAB working solution

### ● BioTnA In Situ RT-PCR Kit with DAB Brown for DIG labeled probe

#### Description

The BioTnA In Situ RT-PCR Kit with DAB is designed to be used for the detection of digoxigenin (DIG)-labeled probes hybridized specific RNA sequence through PCR application in either formalin-fixed, paraffin-embedded tissue or cell samples by chromogenic in situ hybridization (CISH). Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

#### Cat No

TAISP02D-50

TAISP02D-100

#### Pack size

50 reactions

100 reactions

#### Contents :

##### Pre-PCR treatment

Heat Pre-treatment Solution,20x  
Hydrogen Peroxidase Block  
Lysis buffer  
Enzyme Solution

##### PCR reagent

PCR reaction mixture  
Control primer set

##### Post-PCR treatment

Fixation solution  
SSC buffer,20x  
Wash buffer,50X

##### Detection kit

Hybridization buffer  
Sealing gel  
Universal probe DIG labeled  
Control probe DIG labeled  
Hi-Effect Block reagent  
Anti-DIG antibody  
Polymer HRP labeling  
DAB buffer  
DAB chromogen,20x  
Nuclear Blue Solution

##### Slide

Negative control slide\*2  
Control slide\*4


### I. primer design

### II. Sample preparing


- 1.dewax
- 2.rehydration

### III. pretreatment

- 1.De-peroxidase
- 2.Heat per-treatment
- 3.Lysis per-treatment
- 4.Enzyme pre-treatment



### IV. PCR Step



### V. Post-PCR Fixation

### VI. In Situ Hybridization

### VII. Detection

- 1.anti-Dig
- 2.HRP Labeling
- 3.DAB working solution



## ● TUNEL Apoptosis Assay Kit

### Description

The kit can detect fragmented DNA in the nucleus during apoptosis. In this modified TUNEL assay kit, Digoxigenin-nucleotide is labeled at the DNA 3'-OH ends using the natural or recombinant terminal deoxynucleotidyl transferase (TdT or rTdT). Then, mouse anti- DIG and HRP labeled polymer bind to these Digoxigenin nucleotides, which are detected using the peroxidase substrate, 3,3'-diaminobenzidine (DAB), a stable chromogen. Using this procedure, apoptotic nuclei are stained dark brown.

Cat No	Pack size
TAAP01D-20	20 reactions
TAAP01D-50	50 reactions
TAAP01D -100	100 reactions

### Contents :

- Proteinase K
- Permeabilization buffer
- TdT Reaction Buffer
- TdT Enzyme Reagent
- TdT Label Reagent
- Background reducing buffer
- anti- DIG
- HRP labeled polymer
- DAB Chromogen (20x)
- DAB Buffer

## ● TUNEL Apoptosis Assay Kit(FAM)

### Description

The kit can detect fragmented DNA in the nucleus during apoptosis. In this modified TUNEL assay kit, Digoxigenin-nucleotide is labeled at the DNA 3'-OH ends using the natural or recombinant terminal deoxynucleotidyl transferase (TdT or rTdT). The mouse anti- DIG bind to these Digoxigenin nucleotides, which are detected by following the anti-Mouse IgG (H+L)-FAM. The nuclei of apoptotic cells should be observed green color under fluorescence microscope.

Cat No	Pack size
TAAP01F-20	20 reactions
TAAP01F-50	50 reactions
TAAP01F-100	100 reactions

### Contents :

- Proteinase K
- Permeabilization buffer
- TdT Reaction Buffer
- TdT Enzyme Reagent
- TdT Label Reagent
- Background reducing buffer
- anti- DIG
- anti-mouse IgG (H+L)-FAM(488)
- DAPI solution
- Aqua mounting medium

● **Biotna – Tools of molecular Fluorescence multiple stain kit (Provide customized services)**

Cat No	Pack size
TATS01F	1ml

Storage condition : Store at 2-8 °C

Expired date : 18-24 months

**TATS02F**

**Contents :**

Fluorescence Blocking Reagent  
Goat anti-Rabbit IgG (H+L)-488  
Goat anti-Mouse IgG (H+L)-594  
Goat anti-Rabbit IgG-iFluor-680

**TATS03F**

**Contents :**

Fluorescence Blocking Reagent  
Goat anti-Rabbit IgG (H+L)-488  
Goat anti-Rabbit IgG (H+L)-594  
Goat anti-Rabbit IgG-iFluor-680

● **Fluorescence Blocking Reagent**

Cat No	Pack size
TA00B2-15	15ml
TA00B2-50	50ml
TA00B2-100	100ml

● **Goat anti-Mouse IgG (H+L)-594**

Cat No	Pack size
TAFB01-T	500ul

● **Goat anti-Rabbit IgG (H+L)-594**

Cat No	Pack size
TAAB01	500ul

● **Donkey anti-Goat IgG (H+L)488**

Cat No	Pack size
TAFB03-T	500ul

● **Goat anti-Mouse IgG (H+L)-488**

Cat No	Pack size
TAFB01-F	500ul

● **Goat anti-Rabbit IgG (H+L)-488**

Cat No	Pack size
TAFB02-F	500ul

● **Donkey anti-Goat IgG (H+L)-488**


Cat No	Pack size
TAFB03-F	500ul

● **Goat anti-Rabbit IgG-iFluor-680**

Cat No	Pack size
TAFB02-C	500ul


# Histochemistry staining Kit

## ■ Masson Trichrome




Cat No	Major target	volumn
TASS01-125	Collagen fiber (blue), muscle (red)	125ml
TASS01-250		250ml

## ■ Periodic Acid-Schiff(PAS)




Cat No	Major target	volumn
TASS02-125H	Carbohydrates, such as glycogen, glycoprotein, mucus, and basement membrane	125ml
TASS02-250H		250ml
TASS02-500H		500ml
TASS02-250G		250ml
TASS02-500G		500ml

## ■ Alcian blue(pH2.5) 、 Alcian blue(pH1.0)



Cat No	Major target	volumn
TASS03-125pH2.5	acidic polysaccharides, such as cartilage, mucin and others	125ml
TASS03-250pH2.5		250ml
TASS03-125pH1.0		125ml
TASS03-250pH1.0		250ml


## ■ Nissl stain



Cat No	Major target	volumn
TASS04-125	Nissl body	125ml
TASS04-250		250ml


# Histochemistry staining Kit

## ■ Von Kossa




Cat No	Major target	volumn
TASS05-500	Calcium	500ml

## ■ Oil red O (frozen section only,concentrate)




Cat No	Major target	volumn
TASS06	Lipid droplet	500ml

## ■ Toluidine blue




Cat No	Major target	volumn
TASS07	metachromatic dye, commonly used to defined mast cell	125ml

## ■ PicroSirius red



Cat No	Major target	volumn
TASS08-125	Collagen fiber	125ml
TASS08-500		500ml


## ■ Alizarin Red S



Cat No	Major target	volumn
TASS09	Calcium	125ml


# Histochemistry staining Kit

## ■ Safranin O




Cat No	Major target	volumn
TASS10	Cartilage, mucin, mast cell	125ml

## ■ Verhoeff-Van Gieson stain(VVG)




Cat No	Major target	volumn
TASS11	Elastin fiber	125ml

## ■ Congo red-dark




Cat No	Major target	volumn
TASS12	Amyloid substances	125ml

## ■ Prussian blue



Cat No	Major target	volumn
TASS15	Ferric ion (3+)	125ml

## ■ Luxol-Fast blue



Cat No	Major target	volumn
TASS17-125	Myeline	250ml
TASS17-500		500ml

# Ancillaries for Immunohistochemistry

Product	Cat No	Pack size
Hematoxylin for HE	TA01MH-1L	1L
Nuclear blue for IHC (modified hematoxylin)	TA01NB-50	50ml
	TA01NB-100	100ml
	TA01NB-200	200ml
	TA01NB-500	500ml
	TA01NB-1L	1L
Eosin	TA01ES-200	200ml
	TA01ES-1L	1L
Hematoxylin and Eosin	TA01HE-1L	1L
Hydrogen Peroxidase Block	TA00B1-50	50ml
	TA00B1-100	100ml
	TA00B1-200	200ml
	TA00B1-500	500ml
	TA00B1-1000	1000ml
Background Reducing buffer (Immunblock)	TA00B3-15	15ml
	TA00B3-50	50ml
	TA00B3-100	100ml
Background Reducing buffer plus	TA00B4-15	15ml
	TA00B4-50	50ml
	TA00B4-100	100ml
Davidson' s Fixative solution	TADF01-200	200ml
	TADF01-500	500ml
	TADF01-1000	1000ml

Product	Cat No	Pack size
20X Citrate Buffer, pH 6.0	TA00H01-100	100ml
	TA00H01-200	200ml
	TA00H01-500	500ml
	TA00H01-1000	1000ml
20X Tris EDTA Buffer, pH9.0	TA00H02-100	100ml
	TA00H02-200	200ml
	TA00H02-500	500ml
	TA00H02-1000	1000ml
Enzyme Digester I (Pepsin)	TA00E01-15	15ml
	TA00E01-50	50ml
	TA00E01-100	100ml
Enzyme Digester II (proteinase K)	TA00E02-15	15ml
	TA00E02-50	50ml
	TA00E02-100	100ml
Enzyme Digester III (Trypsin)	TA00E03-15	15ml
	TA00E03-50	50ml
	TA00E03-100	100ml
Permont Mount	TA00M1	500ml
Aqua Mount	TA00M2-15	15ml
	TA00M2-50	50ml
	TA00M2-100	100ml
1X PBS, (powder)	TABS01-20	20package



# Ancillaries for Immunohistochemistry

Product	Cat No	Pack size	Product	Cat No	Pack size
10X TBST	TABS02-200	200ml	Goat anti-Mouse IgG (H+L)-HRP	TAAB01	1ml
	TABS02-500	500ml			
	TABS02-1000	1000ml			
20X SSC Buffer	TABS03-200	200ml	Goat anti-Rabbit IgG (H+L)-HRP	TAAB02	1ml
	TABS03-500	500ml			
	TABS03-1000	1000ml			
Antibody Diluent	TABS04-50	50ml	Donkey anti-Goat IgG (H+L)-HRP	TAAB03	1ml
	TABS04-100	100ml			
	TABS04-200	200ml			
	TABS04-500	500ml			
	TABS04-1000	1000ml			
Nuclear Fast Red	TA1NR-50	50ml	Goat anti-Mouse IgG (H+L)-HRP	TAFB02-T	1ml
	TA1NR-100	100ml			
Fast green	TA01FG-100	100ml	Normal Mouse IgG	TANCG01	50ml
	TA01FG-200	200ml			
Methyl Green	TA01MG-100	100ml	Normal Rabbit IgG	TANCG02	50ml
	TA01MG-200	200ml			
Goat probe	TA00C1-15	15ml	Normal Goat IgG	TANCG04	50ml
	TA00C1-50	50ml			
Rat probe	TA00C2-15	15ml	Hydrogen Peroxidase Block	TA00B1-50	50ml
	TA00C2-50	50ml		TA00B1-100	100ml
		TA00B1-200		200ml	
		TA00B1-500		500ml	
		TA00B1-1000		1000ml	
		Fluorescence Blocking Reagent	TA00B2-15	15ml	
			TA00B2-50	50ml	
			TA00B2-100	100ml	
		Background Reducing buffer	TA00B3-15	15ml	
			TA00B3-50	50ml	
			TA00B3-100	100ml	
		Background Reducing buffer plus	TA00B4-15	15ml	
			TA00B4-50	50ml	
			TA00B4-100	100ml	

## ● TnATaq DNA polymerase

### Description

TnATaq DNA polymerase is a thermostable enzyme isolated from *E. coli* which encodes Taq DNA polymerase gene. This enzyme contains 5' -3' polymerase and 5' -3' exonuclease activity.

### Storage buffer

50mM Tris-HCl pH7.9, 50mM KCl, 0.1mM EDTA, 1mM DTT, 0.5mM PMSF, 50% glycerol.

### 10X reaction buffer

Containing 15mM MgCl<sub>2</sub>

### Unit description

One unit is defined as the amount of enzyme that will incorporate 10nmole of dNTP into acid-insoluble material in 30 minutes at 74°C. The reaction conditions are: 50mM Tris-HCl pH8.8, 50mM NaCl, 5mM MgCl<sub>2</sub>, 200uM each of dATP, dCTP, dGTP, dTTP, 10 mg activated calf thymus DNA and 0.1mg/ml BSA in a final volume of 50 ml.

### Storage buffer

50% glycerol (v/v), 20 mM Tris-HCl pH 8.7 at -20°C, 100 mM KCl, 0.1 mM EDTA.

### Source

*E. coli* clone

### Quality control

The enzyme is free of nicking and priming activities, exonucleases and non-specific endonucleases. SDS/PAGE - 95 kD band. Activity and stability tested via thermo-cycling. The error rate per nucleotide per cycle is  $\sim 2.5 \times 10^{-5}$ ; the accuracy is  $\sim 4 \times 10^{-4}$ . Estimated half life at 95°C is 0.5 hours.

### Cat No

TAMB01Z-500  
TAMB01Z-1000  
TAMB01Z-2500

### Pack size

500U  
1000U  
2500U

### Conc.

5U/uL  
5U/uL  
5U/uL

### PCR cycles program


Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	1-3 mins	1
Denaturation	94-95°C	10-60sec	25-35
Annealing	50-68°C	10-30sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

### IMPORTANT:

Annealing temperature should be 2-6°C lower than the primer melting temperature.

### Shipping and Storage conditions

Shipping and temporary storage at -20 and for up to 1 month at room temperature has no detrimental effects on the quality of TnATaq DNA polymerase.





## ● Hot start Taq DNA polymerase

### Description

Hot start Taq DNA Polymerase is designed for Real-Time PCR and Hot- start PCR. It is modified with a special inhibition of PCR at room temperature. This will prevent primer dimers and other artifacts.

### Storage buffer

50mM Tris-HCl pH7.9, 50mM KCl, 0.1mM EDTA, 1mM DTT, 0.5mM PMSF, 50% glycerol.

### 10X reaction buffer

Buffer containing 25mM MgCl<sub>2</sub>

### Unit description

One unit is defined as the amount of enzyme that will incorporate 10nmole of dNTP into acid-insoluble material in 30 minutes at 74°C. The reaction conditions are: 50mM Tris-HCl pH8.8, 50mM NaCl, 5mM MgCl<sub>2</sub>, 200uM each of dATP, dCTP, dGTP, dTTP, 10 mg activated calf thymus DNA and 0.1mg/ml BSA in a final volume of 50 ml.

### Source

E coli clone

### Applications

- Hot Start and real time PCR
- Multiplex PCR
- Amplification of complex genomic and cDNA

### Cat No

TAMB02Z-500  
TAMB02Z-1000  
TAMB02Z-2500

### Pack size

500U  
1000U  
2500U

### Conc.

5U/uL  
5U/uL  
5U/uL

### PCR cycles program


Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	10 mins	1
Denaturation	94-95°C	10-60sec	25-35
Annealing	50-68°C	10-30sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

### IMPORTANT:

Annealing temperature should be 2-6°C lower than the primer melting temperature.

### Shipping and Storage conditions

Shipping and temporary storage at -20 and for up to 1 month at room temperature has no detrimental effects on the quality of Hot start DNA polymerase.



M:100bp ladder

1. Hot start Taq

2. TnATaq

## ● Hi Fidelity polymerase

### Description

Hi Fidelity polymerase is a mixture of thermostable enzymes. It is specifically developed to synthesize length of PCR product up to 25 kb and with low error rate. Hi Fi polymerase synthesizes higher yields of product from genomic DNA, cDNA, and bacterial cultures. It has 2.5 hours half life at 96°C and easily amplify PCR product of G-C rich or secondary structure DNA by adding G- C rich buffer.

### Storage condition

-20°C

### 10X reaction buffer

Buffer containing 25mM MgCl<sub>2</sub>

### Unit description

One unit is defined as the amount of enzyme that will incorporate 10nmole of dNTP into acid-insoluble material in 30 minutes at 74°C. The reaction conditions are: 50mM Tris-HCl pH8.8, 50mM NaCl, 5mM MgCl<sub>2</sub>, 200uM each of dATP, dCTP, dGTP, dTTP, 10 mg activated calf thymus DNA and 0.1mg/ml BSA in a final volume of 50 ml.

### Template

Hi Fidelity Polymerase is suitable for amplifying targets up to 25 kb from the following templates:  
 Genomic DNA: 10–200 ng  
 Plasmid DNA : 1–5 ng  
 cDNA : ~100 ng starting total RNA  
 Amplification of longer targets (up to 15 kb) may be possible, but may require more template and longer elongation times.

### Cat No

TAMB03Z-500  
TAMB03Z-2500

### Pack size

500U  
2500U

### Conc.

5U/uL  
5U/uL

### PCR cycles program


Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	1-3 mins	1
Denaturation	94-95°C	10-60sec	25-35
Annealing	50-68°C	10-30sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

### IMPORTANT:

Annealing temperature should be 2-6°C lower than the primer melting temperature.

### Shipping and Storage conditions

Shipping and temporary storage at -20 and for up to 1 month at room temperature has no detrimental effects on the quality of TnAaq DNA polymerase.



Hi Fidelity polymerase amplified λ DNA from 2.5k to 30k DNA fragment.

## ● Pfu DNA polymerase

### Description

Pfu DNA polymerase is a thermostable enzyme isolated from *Pyrococcus furiosus*. The enzyme replicates DNA at 75°C, catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction. Pfu DNA polymerase possesses 3'→5' exonuclease (proofreading) activity. Base misincorporation is rapidly excised by the proofreading activity of the polymerase. Pfu DNA polymerase is recommended for PCR and primer extension reactions that require high-fidelity. The fragments of Pfu DNA polymerase generated are blunt-ended.

### Error rate

2 x10<sup>-6</sup>

### Unit description

One unit of Pfu DNA Polymerase incorporates 10 nmol of dNTP into acid-insoluble material in 30 min at 74°C.

### Storage buffer

50% glycerol (v/v), 20 mM Tris-HCl pH 8.7 at -20°C, 100 mM KCl, 0.1 mM EDTA.

### Source

E coli clone

### Quality control

The enzyme is free of nicking and priming activities, exonucleases and non-specific endonucleases. SDS/PAGE - 95 kD band. Activity and stability tested via thermo-cycling.

### Cat No

TAMB04Z-500  
TAMB04Z-2500

### Pack size

500U  
2500U

### Conc.

5U/μL  
5U/μL

### PCR cycles program


Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	1-3 mins	1
Denaturation	94-95°C	0.2-1	25-35
Annealing	50-68°C	0.2-1	
Extension	68-75°C	2min/1kb	
Final extension	68-75°C	5	1

### IMPORTANT:

Annealing temperature should be 2-6°C lower than the primer melting temperature.

### Shipping and Storage conditions

Shipping and temporary storage at -20 and for up to 1 month at room temperature has no detrimental effects on the quality of Pfu DNA polymerase.




Pfu DNA polymerase amplified 1kb to 7.5kb DNA Fragment.

## ● 2x PCR mix

### Description

2x PCR mix is optimized mixture. It contains Taq polymerase, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x PCR mix is designed to allow the user to quickly and easily prepare the mixture of reaction. The PCR mix may amplify products up to 3 kb and the products can be cloned into T-vector directly.

Cat No	Pack size	Conc.
TAMB05Z-500	3x1.7 mL	2x
TAMB05Z-2500	5x3x1.7 mL	2x



Comparison of PCR products

## ● 2x Redy mix

### Description

2X Redy Mix is optimized mixture. It contains Taq polymerase, reaction buffer, dNTP, enhancer and red dye as 2-fold concentration. 2x Redy mix is designed to allow the user to quickly and easily prepare the mixture of reaction and ready loading.

Cat No	Pack size	Conc.
TAMB06Z-500	3x1.7 mL	2x
TAMB06Z-2500	5x3x1.7 mL	2x

### Storage condition

long time at -20°C

### Functional Assay

It is tested for performance in the polymerase chain reaction (PCR) to amplify a 500 bp region of the actin cDNA gene. The resulting PCR product is visualized on an ethidium bromide-stained agarose gel.

### PCR cycles program

Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	1-3 mins	1
Denaturation	94-95°C	10-60sec	25-35
Annealing	50-68°C	10-30sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

### Contaminant Assay

No contaminating endonuclease or exonuclease activity detected.

### Application

Amplification of DNA fragments up to 5 kb  
 Labelling of PCR products with modified nucleotides (biotin-dUTP, fluorescein-dUTP)  
 Cycle sequencing  
 Generation of PCR product for TA cloning

### Shipping and Storage conditions

Shipping and temporary storage at -20 and for up to 1 month at room temperature has no detrimental effects on the quality of Taq DNA polymerase.

## ● 2X Hi Fi mix

### Description

Hi Fi polymerase synthesizes higher yields of PCR product from genomic DNA, cDNA, bacterial cultures. It has 2.5 hours half life at 96°C and is easy to amplify PCR product of G-C rich or secondary structure DNA.

2X Hi Fi Mix is optimized mixture. It contains Hi Fi polymerase, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x Hi Fi mix is designed to allow the user to quickly and easily prepare the mixture of reaction. The 2x Hi Fi mix can amplify PCR products up to 10-15 kb and the products can be cloned into T-vector directly.

### 2X Hi Fi Mix

containing 3.6mM MgCl<sub>2</sub>


### Application

Target of DNA fragments is up to 20 kb □Proof reading function  
Amplification of genomic DNA

Cat No	Pack size	Conc.
TAMB07Z-500	3x1.7 mL	2x
TAMB07Z-2500	5x3x1.7 mL	2x

### Storage conditions

long time at -20°C  
short time at 4 °C



## ● 2X Hot start mix

### Description

Hot start Taq DNA Polymerase is designed for Real-Time PCR and Hot-start PCR. It is modified with a special inhibition of PCR at room temperature. This will prevent primer dimers and other artifacts.

2X Hot Start mix is optimized mixture. It contains Hot Start Taq polymerase, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x Hot Start mix is designed to allow the user to quickly and easily prepare the mixture of reaction. The 2x Hot Start mix can amplify PCR products up to 5 kb and the products can be cloned into T-vector directly.

Cat No	Pack size	Conc.
TAMB08Z-500	5x3x1.7 mL	2x
TAMB08Z-2500	5x3x1.7 mL	2x

### Storage conditions

long time at -20°C  
short time at 4 °C

### 2X Hi Fi Mix

containing 3.6mM MgCl<sub>2</sub>


### Application

Target of DNA fragments is up to 3-5 kb  
Hot Start reaction  
Real Time PCR

## ● 2X Pfu mix

### Description

Pfu DNA polymerase is a thermostable enzyme isolated from *Pyrococcus furiosus*. The enzyme replicates DNA at 75°C, catalyzing the polymerization of nucleotides into duplex DNA in the 5'-3' direction. Pfu DNA polymerase possesses 3'-5' exonuclease (proofreading) activity. Base misincorporation is rapidly excised by the proofreading activity of the polymerase. Pfu DNA polymerase is recommended for PCR and primer extension reactions that require high-fidelity. The fragments of Pfu DNA polymerase generated are blunt-ended.



### Cat No

TAMB09Z-500  
TAMB09Z-2500

### Pack size

3x1.7 mL  
5x3x1.7 mL

### Conc.

2x  
2x

### Storage conditions

long time at -20°C

### 2X Pfu Mix

containing 3.6mM MgCl<sub>2</sub>

### Application

Amplification of DNA fragments is up to 7.5 kb Proof reading function.

## ● 5X GC rich buffer

### Description

Improve the amplification of GC-rich targets from genomic DNA or cDNA template. The GC rich buffer can reducing the denature temperature of template.

### Application

GC-rich targets  
Repetitive sequences  
Templates with varying GC content  
Improves fidelity of PCR amplification

### Cat No

TAMB34Z

### Pack size

1 mL

### Storage conditions

long time at -20°C

### Denaturation

The first denaturation step must be more than 5minutes.

## ● 2x qPCR Master mix

### Description

2x qPCR Master mix is designed for quantitative real-time analysis of DNA samples.

- Strong signals and high sensitivity due to fluorescent dye.
- High specificity - no primer dimers, no NTC signal.
- Optimized 2x qPCR Master mixes for different real-time PCR instruments. Master mix formulations are optimized for different machines.

2x qPCR Master mix is ideally suited for:

- Gene expression analysis
- Microarray validation
- Viral load determination qRT-PCR


### Components

2x qPCR Master mix is a 2x mixture of dNTPs, Hotstart Taq polymerase, MgCl<sub>2</sub>, fluorescent detection dye, reference dye (optional), and buffer components.


### Storage condition

long time at -20°C

<b>Cat No</b>	TAMB11Z-500	TAMB11Z-2500
<b>Pack size</b>	3x 1.68 mL 5x3x1.68 mL	3x 1.68 mL 5x3x1.68 mL
<b>Description</b>	2X qPCR master mix with ROX	2x qPCR master mix
<b>qPCR Instruments</b>	ABI,7000,7300,7700,7900, stepOne Plus StepOne™ Eppendorf Realplex 4 ABI7500 Stratagene Mx3000, Mx3005, Mx4000	BioRad CFX96 Roche LightCycler 480 MJ Research Opticon and Opticon 2 MJ Research Chromo 4 Corbett Rotor-gene 600,3000 Eppendorf Realplex 2 Product Application



Comparison: BioTnA and another



The real time PCR results of different concentration of  $\beta$ -actin cDNA

## ● Probe 2X qPCR Master mix

### Description

Probe 2x qPCR Master mix is designed for quantitative real-time analysis of DNA samples.

- Strong signals and high sensitivity due to fluorescent dye.
- High specificity - no primer dimers, no NTC signal.
- Optimized 2x qPCR Master mixes for different real-time PCR instruments. Master mix formulations are optimized for different machines.

Probe 2x qPCR Master mix is ideally suited for:

- Gene expression analysis
- Microarray validation
- Viral load determination qRT-PCR

### Components

Probe 2x qPCR Master mix is a 2X mixture of dNTPs, Hot start Taq polymerase, MgCl<sub>2</sub>, fluorescent detection dye, reference dye (optional), and buffer components.

### Storage condition

long time at -20°C

<b>Cat No</b>	TAMB13Z-500	TAMB13Z-2500
<b>Pack size</b>	3x 1.68 mL 5x3x1.68 mL	3x 1.68 mL 5x3x1.68 mL
<b>Description</b>	Probe 2X qPCR master mix with ROX	Probe 2x qPCR master mix
<b>qPCR Instruments</b>	ABI,7000,7300,7700, 7900, stepOne Plus StepOne™ Eppendorf Realplex 4 ABI7500 Stratagene Mx3000, Mx3005, Mx4000	BioRad CFX96 Roche LightCycler 480 MJ Research Opticon and Opticon 2 MJ Research Chromo 4 Corbett Rotor-gene 600 ,3000 Eppendorf Realplex 2 Product Application



## ● M-MLV reverse transcriptase H-

### Description

MMLV Reverse Transcriptase, from Murine Leukemia virus, is an RNA-dependent DNA polymerase. It synthesizes the cDNA first strand from a single-stranded RNA template to which a primer has been hybridized. MMLV RT also can extend primers hybridized to single-stranded DNA. Second strand cDNA synthesis can be achieved from some RNA templates without an additional DNA polymerase. M-MLV RT (H-) can synthesize 9.5kb products, the largest RNA component in the reaction. However, M-MLV RT synthesized more Full-length cDNA regardless of size.

### Storage condition

long time at -20°C

### Unit description

One unit of activity is the amount of enzyme required to incorporate 1 nmole of dNTP into an acid-insoluble form in 10 minutes at 37°C using polyA-oligo (dT) as template and primer.

### Supplied 5xRT buffer

TrisHCl pH 8.3, KCl, MgCl<sub>2</sub>, DTT

### Storage buffer

50mM Tris-HCl pH8.3, 100mM NaCl, 1mM EDTA, 0.1mM DTT, 0.1% Triton X-100, 50% glycerol

### Cat No

TAMB14Z-100  
TAMB14Z-500

### Pack size

10000 U  
50000 U

### Conc.

200U/μL  
200U/μL

### Activity Assay


200 units of enzyme are used to produce cDNA from 1μg of a 1.2kb control RNA. The minimum specification is 120ng of first strand cDNA made from 1μg of RNA. The cDNA product must be >90% full length.

### Contaminant Activity

To test for endonuclease activity, 1μg of Type I supercoiled plasmid DNA is incubated with 500 units of M-MLV Reverse Transcriptase in 1X Reaction Buffer for one hour at 37° C. Following incubation, the supercoiled DNA is visualized on an ethidium bromide-stained agarose gel to verify the absence of visible nicking or cutting (analysis on 0.4μg of DNA).

### Physical Purity

The purity is >90% as judged by SDS-polyacrylamide gels with Coomassie® blue staining.



Two step RT-PCR

## ● 2X One-tube RT-PCR mix

### Description

The 2x One-tube RT-PCR Hot Start kit is designed for combining two reactions of reverse transcription and PCR. It provides more simple and effective operation in RT-PCR. The reverse transcription step is working by M-MLV RT (H-). The PCR reaction is working by Hot start Taq DNA Polymerase. An especially reaction buffer is provided both for Reverse Transcriptase and Hot Start Taq DNA Polymerase.

### Storage condition

long time at -20°C

### Mix component

1. M-MLV RTase
2. Hot start Taq
3. RT-PCR reaction buffer
4. dNTP
5. stabilizer

### Sensitivity

Targets can generally be detected from < 1 pg to 50 ng polyA RNA(mRNA) or 10 pg to 1 µg total RNA. Even lower amounts of RNA may be successfully amplified by using highly expressed transcripts.

### Cat No

TAMB17Z-50  
TAMB17Z-100  
TAMB17Z-500

### Pack size

50 rx  
100 rx  
500 rx

### PCR cycles program

Step	Temperature	Time	Cycle
RT reaction	37-50°C	30-120mins	1
Initial denaturation	94°C	10mins	1
Denaturation	94°C	0.2-1sec	30-45
Annealing	50-68°C	0.2-1sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

### Analyzing the RT-PCR Products

Analyze the RT-PCR products by 1.0% (w/v) agarose gel electrophoresis. The products will be visible by UV transillumination of the ethidium bromide-stained agarose gel. Store the reaction products at -20°C until needed.

## ● 2X RT Master Mix

### Description

2X RT master mix is designed for convenient and efficient cDNA synthesis. The 2X pre-mixed reagent containing M-MLV RTase, Random 6mers, oligo dT, dNTP mixture and reaction buffer. The reverse transcription reaction just add RNA and water. The protocol is easy and the reaction can be completed in short time.

### Storage condition

long time at -20°C

### 2X RT Master Mix

Containing M-MLV reverse transcriptase, reaction buffer, dNTP, oligo dT primers, random hexamers and stabilizer.

### Application

cDNA synthesis  
PCR screening  
Real-time PCR

### Cat No

TAMB18Z-500  
TAMB18Z-2500

### Pack size

500 uL  
2500 uL

### Conc.

2X  
2X

### Activity Assay


200 units of enzyme are used to produce cDNA from 1µg of a 1.2kb control RNA. The minimum specification is 120ng of first strand cDNA made from 1µg of RNA. The cDNA product must be >90% full length.

### Contaminant Activity

To test for endonuclease activity, 1µg of Type I supercoiled plasmid DNA is incubated with 500 units of M-MLV Reverse Transcriptase in 1X Reaction Buffer for one hour at 37°C.

### Physical Purity

The purity is >90% as judged by SDS-polyacrylamide gels with Coomassie® blue staining.



Lane 1-3 is soybean gene  $\beta$ -actine, Lane 4-6 is E.coli. gene

1: M-MLV RTase + oligo dT

2: 2X master mix

3: M-MLV RTase + specific gene primer

M:100 bp DNA ladder

4: M-MLV RTase + random 6mer

5: 2X master mix

6: M-MLV RTase + specific gene primer

## ● RNAsave solution

### Description

RNAsave solution is a special sample storage liquid, as long as the sample immediately after the fresh sample was immersed in this liquid reagent, RNAsave solution can quickly penetrate into the tissue or other biological sample, the stability and protection of RNA integrity without being degraded, to ensure that the downstream analysis of the dataexpression of the real reaction.

### Cat No

TAMB19Z-500  
TAMB19Z-2500

### Pack size

100 mL  
500 mL

### Storage condition

37°C stable for one day  
18-25°C for 7days  
4 °C for one month  
-20°C for long times

## ● RNase Eraser

### Description

RNase Eraser was formulated using ingredients known to be active against RNase. It effectively removes high levels of RNase contamination. Additionally, Rnase Eraser has been formulated so that it can be used to remove RNase contamination from reaction vessels. The solution is ready to use. If there is a precipitate (as may happen at low temperatures), shake and/or heat at 37°C to bring the precipitate back into solution. The product should be stored at room temperature. Always wear gloves while using RNase Erase™ as prolonged contact with skin may cause irritation. RNase Erase™ is not compatible with corrodible metal surfaces.

### Cat No

TAMB20Z-500

### Pack size

500 mL

### Storage condition

Room temperature

### Instructions

Cleaning work surface  
Cleaning lab apparatus  
Cleaning plastic and glass vessels  
Cleaning pipettes

### Warning

If contact with eyes occurs, immediately flush with water and call a physician. If ingested, do not induce vomiting, give plenty of water and contact a physician.

## ● 5X RNA Protector

### Description

RNA Protector is a mixture of non-toxic reagents for storage and decontamination of RNase for RNA purification. After simple heating step, the contaminated RNase in RNA solution will be eliminated. The treated RNA can be used in all kinds of applications, including northern blot, cDNA synthesis and RT-PCR.

### Cat No

TAMB20Z-500

### Pack size

500 mL

### Storage condition

long time at -20°C

### Usage

Add RNA Protector to final concentration 1X to RNA solution and heat to 60oC for 10-20 minutes.

## ● 100 bp DNA Ladder

### Description

100bp DNA Ladder is suitable for size of linear double stranded DNA fragments from 100 bp to 3 kb. The 13 bands of the ladder contain fragments following Sizes : 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 1500, 2000 bp and 3000 bp. The intensity of 500 bp is about three times than the other bands to serve as reference points.


### Storage buffer

10mM Tris-HCl (pH 8.0), 1mM EDTA

### Storage condition

Long time at -20°C  
At least 6 months at 4°C

Cat No	Pack size	Conc.
TAMB23Z-50	50 ug	0.1 ug/uL
TAMB23Z-250	5x50 ug	0.1 ug/uL



## ● 100 Kb DNA Ladder

### Description

1kb DNA Ladder is suitable for size of linear double stranded DNA fragments from 0.5 kb to 10 kb. The 10 bands of the ladder contain fragments following sizes: 0.25, 0.5, 0.75, 1, 1.5, 2, 3, 4, 5, 6, 8 kb and 10 kb. The intensity of 1 and 4 kb is about three times than the other bands to serve as reference points.


### Storage buffer

10mM Tris-HCl (pH 8.0), 1mM EDTA

### Storage condition

Long time at -20°C  
At least 6 months at 4°C

Cat No	Pack size	Conc.
TAMB24Z-50	100 ug	0.25 ug/uL
TAMB24Z-250	5x100 ug	0.25 ug/uL



0.8% agarose gel analysis

## ● T4 DNA ligase

### Description

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between 5'-phosphate and 3'-hydroxyl ends in double strand DNA or RNA with blunt or cohesive-end termini

### Storage condition

long time at -20°C

### 10X Ligation Buffer

Contain with Tris-HCl, MgCl<sub>2</sub>, DTT, ATP

### Storage buffer

20 mM Tris-HCl (pH 7.5), 1mM DTT, 50mM KCl, 0.1mM EDTA and 50% glycerol

### Inactivation

By heating at 65°C for 10min

### Unit definition

0.01 Weiss unit is defined as the amount of enzyme required to catalyze the ligation of greater than 95% of the Hind III fragments of 1mg of Lambda DNA at 16°C in 20 minutes

### Applications

Joining double-stranded DNA with cohesive or blunt termini and 3' T-A overhangs.

### Cat No

TAMB25Z-50  
TAMB25Z-250

### Pack size


500 U  
2500 U

### Conc.

5 U/μL  
5 U/μL

### Additional information

If the reaction is blunt end or T-A ligation must have more units of T4 DNA ligase. T4 DNA ligase is strongly inhibited by NaCl or KCl if the concentration exceeds 200mM.



Ligation of 1 μg marker using 5U ligase at 25°C

## ● dNTP mix

### Description

The deoxynucleotides (PCR grade) are suitable for many applications where high quality reagents are required. Such procedures include RT, PCR, RT-PCR, DNA labeling reaction and sequencing/ cycle Sequencing analysis.

### Functional assays

Amplification of 8 kb PCR fragment from genomic DNA with Taq DNA polymerase – passed  
 Amplification of 0.6 kb PCR fragment from genomic DNA with Pfu DNA polymerase – passed

#### Cat No

TAMB28Z-S50  
 TAMB28Z-S100  
 TAMB28Z-S250

#### Pack size

4x250 mL  
 4x250 mL  
 4x1 mL

#### Conc.

25mM  
 50mM  
 100mM

### Storage condition

long time at -20°C

### Purity assay

HPLC analysis >99%  
 NMR analysis (inorganic phosphates) – passed  
 Exo-endo deoxyribonucleases test – passed  
 UV analysis – passed  
 Spectrophotometry – passed

## ● dNTP set

### Description

The deoxynucleotides (PCR grade) are suitable for many application where high quality reagents are required. Such procedures include RT, PCR, RT-PCR, DNA labeling reaction and sequencing/ cycle Sequencing analysis.

### Functional assays

Amplification of 8 kb PCR fragment from genomic DNA with Taq DNA polymerase – passed  
 Amplification of 0.6 kb PCR fragment from genomic DNA with Pfu DNA polymerase – passed

#### Cat No

TAMB28Z-25  
 TAMB28Z-100  
 TAMB28Z-250

#### Pack size

4x250 mL  
 4x500 mL  
 4x1000 mL

#### Conc.

25mM  
 50mM  
 100mM

### Storage condition

long time at -20°C

### Purity assay

HPLC analysis >99%  
 NMR analysis (inorganic phosphates) – passed  
 Exo-endo deoxyribonucleases test – passed  
 UV analysis – passed  
 Spectrophotometry – passed

## ● RNA reagent

### Description

RNA REAGENT is a ready-to-use reagent for the isolation of total RNA from cells and tissues. Just a single-step is used during cell homogenization or lysate. 1ml of reagent will treat small quantities of tissue (50-100 mg) or cells (5X 10<sup>6</sup>) well.

### Storage condition

2 to 8°C stable for 12 months

### Application


RT-PCR  
Northern hybridization  
RNase protection d. poly-A+ RNA selection  
Differential display and microarray assay

### Cat No

TAMB29Z-100  
TAMB29Z-200

### Pack size

100 mL  
200 mL



Total RNA analysis

A variety of RNA species has an A260/A280 ratio 1.8 to 2.1 when diluted into TE.

## ● Genomic DNA reagent

### Description

Genomic DNA reagent is a DNA isolation reagent which contains guanidine and detergent mixture. It is a ready to use reagent for the isolation of genomic DNA from various biological sources. During the isolation, a biological sample is lysed (or homogenized) in Genomic DNA reagent and the genomic DNA is precipitated from the lysate with ethanol. Following ethanol wash, DNA is solubilized in water or 8mM NaOH. The procedure can be completed in 10-30 minutes with a 70-100% recovery of genomic DNA.

### Cat No

TAMB30Z-S100  
TAMB30Z-S250

### Pack size

100 mL  
5x100 mL

### Conc.


0.25 ug/uL  
0.25 ug/uL

### Storage condition

Genomic DNA reagent is stable at room temperature for at least two years after the date of purchase.

### Purity assay

Southern analysis  
Dot blot hybridization  
Molecular cloning  
RFLP PCR




1.5 µg DNA analysis



## ● Blue Stain

### Description

Blue Stain is an alternative to traditional Coomassie Blue staining procedures. This ready-to-use stain does not contain methanol, acetic acid or other hazardous solvent. The protein bands are visible directly on gel during the staining process in 10-15 min., then, the gel washing with water yields clear background. The sensitive of detection is up to 10-25ng under standard procedure.



Cat No

TAMB31Z

Pack size

1000 mL

### Storage condition

Room temperature

### Application

Fast: 10-15 can see the protein bands  
Safety: no methanol and acetic acid

## ● Bradford reagent

### Description

The Bradford Reagent is a quick and ready-to-use coomassie-binding colorimetric assay to quantify the total amount of protein. When coomassie dye binds protein in an acidic medium, an immediate shift in absorption maximum occurs from 465 nm to 595 nm with a color change from brown to blue. Note: The concentration of protein should be estimated by reference to absorbance measurements obtained for a series of standard protein dilutions, which are assayed alongside the unknown samples.

Cat No


TAMB32Z

Pack size

500 mL

### Storage condition

The product is stored at 2-8 °C and in an unopened



## ● cDNA Synthesis Ki

### Description

M-MLV cDNA Synthesis kit provides a sensitive and easy-to-use solution for two-step RT-PCR. This kit includes five tubes comprehensive of the reagents required for successful RT-PCR. The M-MLV H-reverse transcriptase is optimized for reliable cDNA synthesis over a wide dynamic range of input RNA. The enzyme is exceptionally well with a wide variety of targets.

### Storage condition

long time at -20°C.

### Kit contents

1. M-MLV RTase H-
2. 5x RT buffer (contain dNTP mix)
3. 5x RNA protector
4. oligodT (15mer)
5. Random primer (6mer)

### Application

1. First-strand cDNA synthesis for subsequent PCR or real-time PCR.
2. RT-PCR validation of gene expression data obtained from microarray experiments.
3. RT-PCR validation and quantification of gene silencing by RNA interference.

### Cat No

TAMB35Z-100  
TAMB35Z-500

### Pack size

100 reaction  
500 reaction

### Activity Assay


200 units of enzyme are used to produce cDNA from 1µg of a 1.2kb control RNA. The minimum specification is 120ng of first strand cDNA made from 1µg of RNA. The cDNA product must be >90% full length.

### Contaminant Activity

To test for endonuclease activity, 1µg of Type I supercoiled plasmid DNA is incubated with 500 units of M-MLV Reverse Transcriptase in 1X Reaction Buffer for one hour at 37° C. Following incubation, the supercoiled DNA is visualized on an ethidium bromide-stained agarose gel to verify the absence of visible nicking or cutting (analysis on 0.4µg of DNA).

### Physical Purity

The purity is >90% as judged by SDS-polyacrylamide gels with Coomassie® blue staining.



Two step RT-PCR

## ● 2x qRT-PCR Master mix

### Description

2x qRT-PCR mix is designed for quantitative real-time analysis of DNA samples.


- Strong signals and high sensitivity due to fluorescent dye.
- High specificity - no primer dimers, no NTC signal.
- Optimized 2x qRT-PCR Master mixes for different real-time PCR instruments. Master mix formulations are optimized for different machines.

### Components

2x qRT-PCR Master mix is a 2X mixture of dNTPs, Hot start Taq polymerase, M-MLV RTase, MgCl<sub>2</sub>, fluorescent detection dye, reference dye (optional), and buffer components.

### Storage condition

long time at -20°C



Quantitative human β-actin gene from 100 ng to 10 pg total RNA

Cat No	Pack size	Description	qPCR Instruments
TAMB36Z-P5 TAMB36Z-P25	3x 1.68 mL 5x3x1.68 mL	2X qRT-PCR master mix with ROX	ABI,7000,7300,7700,7900, stepOne Plus StepOne™ Eppendorf Real-plex 4 ABI7500 Stratagene Mx3000, Mx3005, Mx4000
TAMB36Z-5 TAMB36Z-25	3x 1.68 mL 5x3x1.68 mL	2x qRT-PCR master mix	BioRad CFX96 Roche LightCycler 480 MJ Research Opticon and Opticon 2 MJ Research Chromo 4 Corbett Rotor-gene 600,3000 Eppendorf Realplex 2Product Application

## ● Probe 2x qRT-PCR Master mix

### Description

Probe 2x qRT-PCR mix is designed for quantitative real-time analysis of DNA samples.

- Strong signals and high sensitivity due to fluorescent dye.
- High specificity - no primer dimers, no NTC signal.
- Optimized 2x qPCR Master mixes for different real-time PCR instruments. Master mix formulations are optimized for different machines.

### Components

Probe 2x qRT-PCR Master mix is a 2X mixture of dNTPs, Hot start Taq polymerase, M-MLV RTase, MgCl<sub>2</sub>, reference dye (optional), and buffer components.

### Storage condition

long time at -20°C



Cat No	Pack size	Description	qPCR Instruments
TAMB37Z-P5 TAMB37Z-P25	3x 1.68 mL 5x3x1.68 mL	Probe 2X qRT-PCR master mix with ROX	ABI,7000,7300,7700,7900, stepOne Plus StepOne™ Eppendorf Realplex 4 ABI7500 Stratagene Mx3000, Mx3005, Mx4000
TAMB37Z-5 TAMB37Z-25	3x 1.68 mL 5x3x1.68 mL	Probe 2x qRT-PCR master mix	BioRad CFX96 Roche LightCycler 480 MJ Research Opticon and Opticon 2 MJ Research Chromo 4 Corbett Rotor-gene 600,3000 Eppendorf Realplex 2Product Application

## ● HyFect™ DNA Transfection Reagent

### Description

HyFect™ DNA Transfection Reagent is a low cell toxicity and non-liposome reagent optimized for DNA delivery. It provides effectively, reproducibly, and affordable benefits for scientific research.

### Protocol

- 1. Cell preparation:** Cells should be seeded before 16 to 20 hours prior to transfection with around 70% confluency. The medium should be refreshed 30 min before transfection. Usually, culture media with serum does not affect transfection.
- 2. DNA preparation:** DNA plasmid for transfection should be with high purity (A260/A280 = 1.8-1.9) to ensure efficient transfection mixture preparation.
- 3. Mixture preparation:** Guideline of DNA plasmid and HyFect™ DNA Transfection Reagent amount and ratio refer to table 1. In brief, plasmid and transfection reagent are diluted in serum-free culture medium for 5 min then mixed together and tip gently for another 25 min.
- 4. Transfection:** Add mixtures into cell culture dish / plate. The mixture could be removed after 6 to 48 hours and refilled with culture medium.

Cat No

TAMB33L

Pack size

500 µL

### Storage condition

Store at -20 °C for long term

**Table 1. Recommended formula of transfection mixture**

Culture Dish/Plate	Media Volume	Plasmid/ Serum-free medium	HyFect™ DNA Transfection Reagent/ Serum-free medium
96-well	100 µL	250 ng/ 10 µL	0.75 µL/ 10 µL
24-well	500 µL	500 ng/ 25 µL	1.5 µL/ 25 µL
12-well	700 µL	750 ng/ 35 µL	2.25 µL/ 35 µL
6-well	1 mL	1 µg/ 50 µL	3 µL/ 50 µL
6 cm	3 mL	2.5 µg/ 150 µL	7.5 µL/ 150 µL
10 cm	6 mL	5 µg/ 300 µL	15 µL/ 300 µL

## ● DNA Safe Stains

### Description

DNA Safe Stains is a safe nucleic acid stain to replace traditional ethidium bromide (EtBr, a potent mutagen) for detecting double-stranded DNA, single-stranded DNA, and RNA in agarose gel. It emits green fluorescence when bound to dsDNA and red fluorescence when bound to ssDNA or RNA. DNA Safe Stains products are non-carcinogenic by the Ames-test (negative results in both mouse marrow chromophilous erythrocyte micronucleus and mouse spermary spermatocyte chromosomal aberration tests).

### Protocol

1. Prepare hot melting agarose gel solution (concentration from 0.5-2.5%).
2. Add 1  $\mu\text{L}$  of DNA Safe Stains to every 20 mL of gel solution. Swirl gently to mix the solution and avoid forming bubbles.
3. Allow the agarose gel to cool until solidified. Load samples on the gel and perform electrophoresis.
4. Visualize the bands under UV illumination.

Cat No

TAMB22L

Pack size

1 mL

### Storage condition

Store at 4 °C

### Note

1. The thickness of gel should be less than 0.5 cm since thick gels may decrease sensitivity.
2. Repeated melting of gels containing DNA Safe Stains may result in low sensitivity.
3. DNA Safe Stains may irritate skin and eyes. Please wear gloves while handling.



