



Product name: Package 1 Catalog No.: TAPS01

Introduction:

The scanning parameters for the above staining were adjusted based on slides prepared using BioTnA reagents and staining conditions.

Kit Contents:

Kit Contents	Advanced Image Settings	
HE_Rat-TMA		
	Gamma correction Gamma factor: 0.917 Pixel intensity deviation: 10	Image enhancement Minimum input: 8 Maximum input: 239
Exposure: 139µm Image adjustment Brightness: 50 Contrast: 0	Color adjustment Red: 1.008 Green: 1.008 Blue: 1.040	Color enhancement Correction factor: 0 Sharpness Sharpness factor: 1
Masson Trichrome Pig Skin	Gamma correction Gamma factor: 0.782 Pixel intensity deviation: 12	Image enhancement Minimum input: 13 Maximum input: 236
Exposure: 145µm Image adjustment Brightness: 50 Contrast: 0	Color adjustment Red: 1.048 Green: 1.040 Blue: 1.000	Color enhancement Correction factor: 1 Sharpness Sharpness factor: 1



Picro Sirius Red Pig-Skin		
	Gamma correction	Image enhancement
	Gamma factor: 0.764	Minimum input: 12
e e till a di de de dille	Pixel intensity deviation: 11	Maximum input: 238
	·	
	Color adjustment	Color enhancement
Exposure: 145µm	Red: 1.008	Correction factor: 1
Image adjustment	Green: 1.040	Sharpness
Brightness: 50	Blue: 1.024	Sharpness factor: 1
Contrast: 0		-
Safranin O_Mouse-Knee		
	Gamma correction	Image enhancement
	Gamma factor: 0.635	Minimum input: 8
	Pixel intensity deviation: 9	Maximum input: 239
	Color adjustment	Color enhancement
Exposure: 145µm	Red: 1.104	Correction factor: 1
Image adjustment	Green: 1.080	Sharpness
Brightness: 50	Blue: 1.000	Sharpness factor: 1
Contrast: 0		
IHC(a-SMA)_Dog-Tumor		
STATE OF BRIDE OF ANGLES OF THE STATE OF THE	Gamma correction	Image enhancement
	Gamma factor: 0.723	Minimum input: 19
	Pixel intensity deviation: 14	Maximum input: 231
	Color adjustment	Color enhancement
Exposure: 145µm	Red: 1.000	Correction factor: 1
Image adjustment	Green: 1.048	Sharpness
Brightness: 50	Blue: 1.080	Sharpness factor: 1
Contrast: 0		