





Package 4 

Product name: Package 4
 Catalog No.: TAPS04

Introduction:

The scanning parameters for the above staining were adjusted based on slides prepared using BioTnA reagents and staining conditions.

Kit Contents:

Kit Contents	Advanced Image Settings
<p>Oil Red O</p>  <p>Exposure: 145µm Image adjustment Brightness: 50 Contrast: 0</p> 	<p>Gamma correction Gamma factor: 0.653 Pixel intensity deviation: 12</p> <p>Image enhancement Minimum input: 13 Maximum input: 226</p>
<p>Alcian Blue_Mouse-Intestine</p>  <p>Exposure: 152µm Image adjustment Brightness: 50 Contrast: 0</p> 	<p>Gamma correction Gamma factor: 0.717 Pixel intensity deviation: 12</p> <p>Image enhancement Minimum input: 6 Maximum input: 227</p>
	<p>Color adjustment Red: 1.000 Green: 1.032 Blue: 1.040</p> <p>Color enhancement Correction factor: 0 Sharpness Sharpness factor: 1</p>
	<p>Color adjustment Red: 1.000 Green: 1.080 Blue: 1.064</p> <p>Color enhancement Correction factor: 1 Sharpness Sharpness factor: 1</p>